

# AFFINITY PURIFICATION OF SOLUBLE HIS-TAGGED PROTEINS

[www.ptglab.com](http://www.ptglab.com)

---

1.

## Lyse cells:

- a. Suspend the cell pellet in 30–35 ml of His-washing buffer with 10 mM PMSF.
  - b. Sonicate cells in an ice-bath at 200 W for 6 min.
  - c. Rotate the lysed solution for 1 h at 4°C.
  - d. Centrifuge the cell lysate for approximately 13 min at 8000 rpm, 4°C.
- 

2.

## Bind protein to beads:

- a. Transfer the supernatant to 600 µl of His-beads.
  - b. Rotate the mixture overnight at 4°C.
  - c. Collect the beads by centrifugation at 2000 rpm for 10–30 seconds, 4°C. The protein-bound beads are collected in eppendorf tubes.
- 

3.

## Wash away the unbound proteins from beads:

- a. Wash the beads 3 times with 1 ml of His-washing buffer. Discard the supernatant.
- 

4.

## Elute proteins from beads:

- a. Add 300 µl of His-elution buffer to the beads.
- b. Rotate the mixture for 1 h at 4°C.
- c. Collect the supernatant by centrifugation at 300 rpm for 10–30 seconds.
- d. Repeat steps 4 a–c.
- e. Combine the eluent (total volume of 600 µl).
- f. Check the molecular weight and purity of the enriched protein by SDS-PAGE analysis.

# AFFINITY PURIFICATION OF SOLUBLE HIS-TAGGED PROTEINS

[www.ptglab.com](http://www.ptglab.com)

## Buffers Needed

<b>His-washing buffer</b>	<b>1000 ml</b>
20 mM Imidazole	1.36 g
1x PBS buffer	1000 ml
Adjust to pH 7.0	
<b>His-elution buffer</b>	<b>1000 ml</b>
300 mM Imidazole	20.42 g
10% Glycerol	100 ml
1x PBST buffer	900 ml
Adjust to pH 7.0	
<b>GST-washing buffer (PBST buffer)</b>	<b>1000 ml</b>
58 mM Na <sub>2</sub> HPO <sub>4</sub>	8.24 g
17 mM NaH <sub>2</sub> PO <sub>4</sub>	2.04 g
68 mM NaCl	3.98 g
Add ddH <sub>2</sub> O to 1000 ml	
Adjust to pH 7.4	