


Nano-BiTEs: Bispecific T cell engagers based on nanobodies

Christian Linke-Winnebeck^{[1]*}, Michael Metterlein^[1], Larisa Yurlova^[1], Sabrina Wendler^[1], Simge Yüz^[2], Björn Tränkle^[3], Jacqueline Bogner^[1], Tanja Ertl^[1], Felix Hartlepp^[1]


^[1]ChromoTek GmbH, Planegg-Martinsried, Germany; ^[2]Tumorbiology, Natural and Medical Sciences Institute, Reutlingen, Germany; ^[3]Pharmaceutical Biotechnology, University of Tübingen, Reutlingen, Germany; *Email to: c.linkewinnebeck@chromotek.com

Based on our long-standing experience in the field of nanobodies (V_HH), we have developed Nano-BiTEs. These bispecific T cell engagers are fusions of a nanobody against the tumor antigens Her2 or CEA and an scFv against the T cell receptor (TCR) subunit CD3ε. Nano-BiTEs are designed as reagents for the efficient activation of T cells in the context of tumour cells.




Nano-BiTE format and affinity

Nano-BiTE constructs



Format of Nano-BiTEs:

- Fusion of V_HH and scFv
- scFv anti-human CD3ε




Production:

- Mammalian expression
- Affinity chromatography
- Endotoxin-free

High affinity of the Nano-BiTEs

Biolayer interferometry (BLI)




	Ligand	K_D (nM)	k_a (M ⁻¹ s ⁻¹)	k_d (s ⁻¹)
VHH0219	α-Her2×α-CD3	3	$2.5 \cdot 10^5$	$7.8 \cdot 10^{-4}$
	Her2	1	$7.2 \cdot 10^5$	$7.3 \cdot 10^{-4}$
VHH0221	α-CEA×α-CD3	12	$4.2 \cdot 10^5$	$5.0 \cdot 10^{-4}$
	CEA	0.8	$5.3 \cdot 10^5$	$4.1 \cdot 10^{-4}$


→ High affinity of Nano-BiTEs will allow targeted tumour and T cell binding.

Antigen-binding of Nano-BiTEs

Concomitant binding to CD3ε and tumour antigen *in vitro* (BLI)




Experimental set-up




Binding to CD3ε and tumour antigens on cells

Immunofluorescence microscopy



→ Nano-BiTEs recognise their antigens CD3ε and Her2 or CEA *in situ* on cells.

Flow cytometry



Conclusion

Our Nano-BiTEs are a novel tool format to activate T cells specifically by recruiting them to tumour cells. Nano-BiTEs can be used in immuno-oncology as standards or to achieve a defined level of T cell activation to quantify the effect of novel drug candidates.

Acknowledgements:
We are grateful to Stefanie Urlinger (iOmx Therapeutics) for fruitful discussions.

References:

[1] Zheng, L. et al. 2019: Structural basis of assembly of the human T cell receptor-CD3 complex. Nature 573, pp. 546-552. PDB 6JXR. [2] Garboczi, D. N. et al. 1996: Structure of the complex between human T-cell receptor, viral peptide and HLA-A2. Nature 384, pp. 134-141. PDB 1AO7. [3] Kjer-Nielsen, L. et al. 2004: Crystal structure of the human T cell receptor CD3ε heterodimer complexed to the therapeutic mAb OKT3. PNAS 101, pp. 7675-7680. PDB 1SY6.

Disclaimer:

ChromoTek is a registered trademark of ChromoTek GmbH. Nanobody is a registered trademark of Ablynx, a Sanofi company. Alexa Fluor is a registered trademark of Life Technologies Corporation, a part of Thermo Fisher Scientific Inc.