

## RBC Lysis Buffer (10X)

1. At room temperature (18-25°C), dilute 1 part of 10× red cell lysis concentrate with 9 parts deionised water to obtain 1× working solution.
2. Following incubation under conditions specified in the antibody manual, proceed with red cell lysis.
3. Gently mix whole blood before red cell lysis.
4. Add 2 mL working solution to each 100 uL fresh anticoagulated blood sample. Invert to mix thoroughly, then lyse at room temperature for 5 minutes.
5. Centrifuge at 300-400 g for 5 minutes at room temperature. Discard the red supernatant. (If incomplete lysis occurs, repeat step 2. Normally, slight amounts of red blood cells do not affect subsequent testing.)
6. Add 2 mL PBS, mix thoroughly, then centrifuge at 300-400 g for 5 minutes at room temperature. Discard the supernatant.
7. Add 200 uL PBS, mix thoroughly, and proceed to instrument analysis promptly.

---

### Proteintech Group, Inc.

*(HQ in North America)*

5500 Pearl Street, Ste 400  
Rosemont, IL 60018, USA  
P +1 (888) 478-4522 (toll free in USA)  
F +1 (773) 598-4321  
+1 (847) 928-2189  
E [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

### Proteintech Europe

P +44 161 839 3007  
F +44 161 241 3103  
E [europe@ptglab.com](mailto:europe@ptglab.com)

### Proteintech China

E [proteintech-CN@ptglab.com](mailto:proteintech-CN@ptglab.com)

### Proteintech Singapore

P +65 9785 0861  
E [singapore@ptglab.com](mailto:singapore@ptglab.com)