

Flow Cytometry Phosphorylated Protein Fixation/Permeabilization Kit

This protocol was developed for analysis of 1×10^6 cells.

1. Bring all buffers to Working Temperature prior to use (see Kit Components Table).
2. Wash cells with PBS. Discard supernatant.
3. Add 0.5 mL Fixation Buffer. Mix well and incubate at RT for 15 min.
4. Centrifuge at 400-600 g for 5 minutes and aspirate (do not pour off) the supernatant.
5. Add 0.5 mL pre-cooled (-20°C) Permeabilization Buffer. Mix well by pipetting and incubate at RT for 5 minutes.
6. Centrifuge at 400-600 g for 5 minutes and aspirate (do not dump) the supernatant.
7. Dilute the Staining Buffer (10X) to 1X with deionized water.
8. Add 1.5 mL Staining Buffer (1X). Mix well and centrifuge at 400-600 g for 5 minutes and aspirate (do not pour off) the supernatant.
9. Repeat wash step 8.
10. Resuspend cells in 100 μL Staining Buffer (1X).

Proteintech Group, Inc.

(HQ in North America)

5500 Pearl Street, Ste 400
Rosemont, IL 60018, USA
P +1 (888) 478-4522 (toll free in USA)
F +1 (773) 598-4321
+1 (847) 928-2189
E proteintech@ptglab.com

Proteintech Europe

P +44 161 839 3007
F +44 161 241 3103
E europe@ptglab.com

Proteintech China

E proteintech-CN@ptglab.com

Proteintech Singapore

P +65 9785 0861
E singapore@ptglab.com

Flow Cytometry Phosphorylated Protein Fixation/Permeabilization Kit

11. Stain cells with desired antibodies at pre-titrated concentrations.
12. After incubation, wash cells twice as in step 8.
13. If using secondary antibody, dilute the secondary antibody using Staining Buffer (1X).
Add 100 uL of diluted secondary antibody and incubate.
14. After incubation, wash cells twice as in step 8.
15. Resuspend cells in 200 uL Staining Buffer (1X).
16. Analyze sample with flow cytometer.

Proteintech Group, Inc.

(HQ in North America)

5500 Pearl Street, Ste 400

Rosemont, IL 60018, USA

P +1 (888) 478-4522 (toll free in USA)

F +1 (773) 598-4321

+1 (847) 928-2189

E proteintech@ptglab.com

Proteintech Europe

P +44 161 839 3007

F +44 161 241 3103

E europe@ptglab.com

Proteintech China

E proteintech-CN@ptglab.com

Proteintech Singapore

P +65 9785 0861

E singapore@ptglab.com