

# Spot-Label ATTO 488

Only for research applications, not for diagnostic or therapeutic use

**1. Introduction** Small peptide tags are useful for the labelling and detection of proteins using immunostaining, immunoblotting, or immunoprecipitation techniques. The ChromoTek Spot-Tag® is a short 12 amino acid affinity tag (sequence: PDRVRAVSHWSS), which can be cloned either N- or C-terminally to a protein of interest. This tag can be efficiently immunostained with the novel Spot-Label affinity reagent. The Spot-Label consists of a small recombinant bivalent alpaca single-domain antibody fragment covalently conjugated to a fluorescent dye. It enables the fluorescence-based Western-blot detection and immunofluorescence microscopy analysis of Spot-Tag fusion proteins. Due to the small size of the Spot-Label, immunostaining of the Spot-Tag with the Spot-Label minimizes the “linkage error” for super-resolution microscopy applications (e.g. STED and dSTORM). In addition, the Spot-Label has a superior tissue penetration rate, better access to the Spot epitope, and higher labelling density.

## 2. Content

Reagent	Quantity	Code
Spot-Label ATTO 488	50 µL	eba488-50
Spot-Label ATTO 488	10 µL	eba488-10

## 3. Properties

### Description:

Recombinant alpaca single-domain antibody for the analysis of Spot-Tag fusion proteins.

### Specificity:

This antibody fragment is reactive against the Spot-Tag (PDRVRAVSHWSS).

### Product Type:

Primary antibody, conjugated to ATTO 488

### Isotype:

V<sub>H</sub>H (Nanobody), alpaca monoclonal, bivalent

### Purity:

Affinity-purified antibody fragment

### Form:

Liquid

### Storage Buffer:

Buffered aqueous solution (PBS)

### Preservative:

0.09% Sodium Azide

### Safety datasheet (SDS) for this product:

**Sodium Azide SDS**

### Concentration:

1 g/L

### Optical properties:

**ATTO 488:** Excitation range 480 - 510 nm ( $\lambda_{abs}$ = 501 nm)

Emission range 520 - 560 nm ( $\lambda_{fl}$ = 523 nm)

For further information please refer to <http://www.atto-tec.com>

#### 4. Stability and Storage

Shipped at ambient temperature. Upon receipt store at +4°C / 40°F.  
Stable for 6 months. Do not freeze. Protect from light.

#### 5. IF Protocol

1. **Fixation:** Fix cells seeded on coverslips in 3.7% formaldehyde in PBS for 10 min at room temperature.  
*Note: Always prepare a fresh formaldehyde dilution.*
2. Wash samples three times with PBS (Phosphate Buffered Saline). Do not store fixed cells.
3. **Permeabilization:** Add PBS containing 0.5% Triton X-100 to samples and incubate for 5 min at room temperature.  
*Note: Alternatively, use ice-cold 100% methanol for permeabilization.*
4. Wash samples twice with PBS.
5. **Blocking:** Add 4% BSA in PBS to samples and incubate for 20 min at room temperature.  
*Note: If necessary, use additional blocking reagents (e.g. 10% normal serum in PBS or Image-iT™ FX Signal Enhancer from ThermoFischer Scientific) and extend the blocking time up to 60 min.*
6. **Spot-Label incubation:** Dilute Spot-Label 1:1,000 - 10,000 in blocking buffer and incubate for overnight at +4°C.  
*Note: For multiplexing protocols, you can combine Spot-Label with another primary or secondary antibody.*
7. Wash samples three times for 5-10 min in PBS.
8. If required, counterstain with DNA fluorescent dyes, e.g. DAPI in PBS. Proceed with imaging directly or mount samples, if necessary.
9. **Mounting:** Rinse sample briefly in water to prevent salt crystal formation. Mount in ProLong™ Diamond Antifade Mountant from ThermoFischer Scientific or other mounting media with anti-fading agents.

#### Suggested buffer composition

Buffer	Composition
Fixation buffer	3.7% formaldehyde in PBS
Permeabilization buffer	PBS; 0.5% Triton X-100
Wash buffer	PBS
Blocking buffer	4% BSA (w/v); PBS

#### 6. Western blot

1. **Preparation:** Separate your sample of interest on an SDS-PAGE gel and transfer onto a nitrocellulose membrane according to standard protocols.
2. **Blocking:** Incubate membrane with 5 % milk powder in PBS or TBS + 0.075 % Tween-20 (PBST or TBST).
3. **Spot-Label incubation:** Dilute Spot-Label in 5 % milk powder in PBST or TBST. The recommended starting dilution is 1:5,000. Add diluted Spot-Label to membrane and incubate at 4 °C overnight.  
*Note: The optimal dilution depends on the application and should be determined by the user. A titration from 1:1,000 up to 1:20,000 is recommended.*

4. Wash samples three times for 5-10 min in PBS or TBST.
5. **Detection:** Image fluorescence using a fluorescence scanner or similar and using appropriate settings

*Note: Preprogrammed imaging settings for fluorescein, AlexaFluor 488, or Cy2 will work also for ATTO 488.*

**Support/  
Troubleshooting**

Please refer to our FAQ section at [www.chromotek.com](http://www.chromotek.com) or contact [support@chromotek.com](mailto:support@chromotek.com)

**Related Products**

Spot-Tag Toolbox	Code
Spot-Trap® Agarose	eta-20
Spot-Trap® Magnetic Agarose	etma-20
Binding control agarose beads	bab-20
Binding control magnetic agarose beads	bmab-20
Spot V <sub>H</sub> H, recombinant binding protein	etb-250
Spot-Label ATTO 594	eba594-50
Spot-Tag peptide	ep-1
Spin columns	sct-10; sct-20; sct-50

**Trademarks and  
Copyrights**

Spot-Label, Spot-Tag and Spot-Trap are registered trademarks of ChromoTek GmbH. *Image-iT*<sup>™</sup> and ProLong<sup>™</sup> are trademarks of Thermo Fischer Scientific Inc.