For Research Use Only

SignalBright Pro Chemiluminescent Substrate



Catalog Number: PK10011

Description

SignalBright Pro is a luminol-based, chemiluminescent substrate for HRP (horseradish peroxidase) which provides a much higher sensitivity over traditional ECL reagents. SignalBright Pro provides mid to high-femtogram detection levels in Western Blot and ELISA analysis, with long lasting signal up to 24 hours. Key Characteristics:

- · Sensitivity: Mid to high femtogram
- · Signal duration: >12 hours
- · Compatibility: suitable for PVDF and nitrocellulose membranes
- · Comparable sensitivity to SuperSignal™ West Dura.

Instructions for Use

SignalBright Pro Chemiluminescent Detection is prepared by combining a 1:1 ratio of the **Chemiluminescent Substrate Solution** and **Stable Peroxide Solution** to create a working solution. Allow the working solution to come to room temperature before applying to the membrane, make sure this is done in darkness. Apply a minimum of 0.1ml of the working solution per membrane cm2. Leave 5 minutes for signal to develop prior to imaging.

- Recommended primary antibody concentration: 0.02-1 ug/ml
- Recommended HRP-conjugated secondary antibody concentration: **4-20** ng/ml

Package

Product Information


20mL/100mL

| Cat. | Size | Content |
|---------|-------|--|
| PK10011 | 20mL | 10mL Chemiluminescent Substrate Solution |
| | | 10mL Stable Peroxide Solution |
| | 100mL | 50mL Chemiluminescent Substrate Solution |
| | | 50mL Stable Peroxide Solution |


Protocol

SignalBright Pro Chemiluminescent Substrate Protocol

Validation Data



Serial dilutions of RUVBL primary antibody (10210-2-AP)
br>Primary: Proteintech RUVBL1 (10210-2-AP) at various dilutions (see image)
br>Secondary: Proteintech HRP-conjugated Affinipure Goat Anti-Rabbit IgG (SA00001-2); 1:6,000
br>Exposure time: 1 min
br>Chemiluminescent substrates from leading competitor, SignalBright Pro (PK10011),...
SignalBright Max (PK10013)



Serial dilutions of HEK-293 cell lysates https://doi.org/10.210-2-AP); 1:8,000 https://doi.org/10.210-2-AP); 1:6,000 https://doi.org/10.210-2-AD); 1:6,000 https://doi.org/10.210-2-AD); 1:6,000 <a href="https://doi.org/10.210-2-AD); 1:6,000 https://doi.org/10.210-2-AD); 1:6,000 <a href="https://doi.org/10.210-2-AD); 1:7,000 <a href="https://doi.org/10.210-2-AD</a