

CoraLite® 488-594 TUNEL Apoptosis

Instructions

Materials required but not supplied with kit:

1. PBS buffer (pH~7.4)
2. 4% paraformaldehyde (in PBS)
3. Bovine serum albumin (BSA) or normal sheep/bovine serum
4. 70% ethanol (optional)
5. Dewaxing solvent (if using paraffin section samples)

Cell samples

1. Prepare a negative control sample (add TUNEL reaction solution without TdT enzyme).
2. Wash the cells twice with PBS.
3. Add 4% paraformaldehyde (pH 7.4) solution in excess and incubate at 4°C for 30 min.
4. Wash the cells twice with PBS.
5. Add 70% ethanol pre-cooled on ice and incubate at -20°C for 4 h. Cells can be stored for one week in 70% ethanol at -20°C. Alternatively, the cells can be permeabilized with 0.2% Triton X-100 in PBS solution and incubated at room temperature for 20 minutes.
6. Wash the cells twice with PBS.

Paraffin tissue sections

- 1. Immerse slides in xylene for 5 minutes. Repeat this step in fresh xylene for a further 5 minutes.**
Note: Xylene is toxic and volatile. Please perform this in a fume hood.
- 2. Immerse slides in absolute ethanol for 5 minutes. Repeat this step for a further 5 minutes in fresh ethanol.**
- 3. Rehydrate sections by sequentially incubating with 100%, 95%, 80%, and 60% ethanol for 5 minutes each.**
- 4. Immerse slides in distilled water for 3 minutes each. Use filter paper to carefully absorb the excess liquid around the sample.**
- 5. Use an immunohistochemical pen to trace around the samples.**

- 6. Dilute Proteinase K solution 1:100 in 1mL of PBS to a final concentration of 20 µg/mL. Add 20 µg/mL Proteinase K solution to each sample in a dropwise manner until the solution covers the entire area of the sample (~100 µL), followed by incubation for 20min at room temperature. (The incubation time and temperature of Proteinase K should be optimized according to your tissue samples.)**
Note: Over-incubation with Proteinase K may cause the sections to fall off, so optimize the length of the incubation time. Typically, the time required is 10-30 minutes. For a section of around 4 µm: 10 minutes, and for sections around 30 µm: 30 minutes.
- 7. Wash the sections with PBS twice for 5 minutes each time, absorb the excess liquid with filter paper, and place the processed sample in a humid box to keep the sample moist.**
Note: Proteinase K must be thoroughly washed off during this step to prevent interference with subsequent labeling.

Frozen tissue sections

1. Place the frozen sections on a rack at room temperature for 20 minutes to thaw.
2. Immerse the slides in 4% paraformaldehyde solution (in PBS) and fix at room temperature for 30 min.
3. Wash the slides twice with PBS for 5 minutes each time.
4. Use filter paper to carefully dry the liquid around the sample on the slide.
5. Dilute Proteinase K solution 1:100 in 1mL of PBS to a final concentration of 20 µg/mL. Add 20 µg/mL Proteinase K solution to each sample in a dropwise manner until the solution covers the entire area of the sample (~100 µL), followed by incubation for 20min at room temperature.

Note: Over-incubation with Proteinase K may cause the sections to fall off, so optimize the length of the incubation time. Typically, the time required is 10-30 minutes. For a section of around 4 µm: 10 minutes, and for sections around 30 µm: 30 minutes.

6. Wash the sections with PBS twice for 5 minutes each time, absorb the excess liquid with filter paper, and place the processed sample in a humid box to keep sample moist.

Note: Proteinase K must be thoroughly washed off during this step to prevent interference with subsequent labeling.

Positive control treatment (*only your positive control sample needs to perform this step, your main experimental samples directly undergo the TUNEL reaction in step 5*)

1. Dilute the 10 X DNase I Buffer with ddH₂O in a ratio of 1:10 to make 1 X DNase I Buffer as a working concentration.
2. Add 100 µL 1 X DNase I Buffer dropwise to your sample, and incubate at room temperature for 5 minutes.
3. Dilute DNase I (2 U/µL) 1:100 with 1 X DNase I Buffer for working solution with a final concentration of 20 U/mL.
4. Gently aspirate the excess liquid. Next add 100 µL of 20 U/mL DNase I working solution to your slides in a dropwise manner, and incubate at room temperature for 10 min.
5. Gently aspirate the excess liquid, and wash the sample twice with PBS.

TUNEL reaction

- 1. Prepare a master mix of the TUNEL reaction mixture in advance:** for each sample allow 50 μ L TUNEL reaction buffer with 1 μ L TdT enzyme.
- 2. Add 100 μ L of the equilibration buffer to each sample and incubate for 5 minutes.**
- 3. Discard the equilibration buffer, carefully aspirate the excess liquid around the sliced samples with filter paper, and add 50 μ L of the TUNEL reaction mixture to each sample.**
 - For adherent cells, evenly cover the sample with the buffer solution before carefully placing a glass coverslip. Incubate the samples in the dark at 37°C for 60 min.*
 - Suspended cells can be added to a microtiter plate and incubated on a plate shaker, alternatively you can gently shake the reaction tube every 15 minutes. Incubate in the dark at 37°C for 60 min.*
 - For tissue samples, cover the sample evenly with buffer and carefully place a glass cover slip. Place the sample in a humid box and incubate at 37°C for 2 hours. Spread a paper towel with a small amount of water on the bottom of the wet box to maintain humidity. Incubate for 2 h at 37°C in the dark.*
- 4. Gently aspirate the reaction solution. Next, soak and rinse twice for 5 mins in a 1X PBS staining tank. Then use an appropriate**

amount of 0.1% Triton X-100 in PBS containing 5 mg/mL BSA buffer to wash the sample 3 times, 5 min each time to reduce background.

- 5. Counterstaining (optional):** Add 2 μ g/mL DAPI dropwise to each sample and incubate for 10 min at room temperature in the dark. Next, remove stain and rinse the slides 3 times in 1X PBS for 5 minutes each time.
- 6. Mounting (optional):**
 - Immerse the slides in distilled water for 5 minutes.*
 - Dehydrate sections by sequentially incubating with 100%, 95%, 80%, and 60% ethanol for 5 minutes each.*
 - Then immerse slides in fresh xylene for 5 minutes and repeat this step for another 5 minutes in fresh xylene. (Note: xylene fumes are toxic so perform this step in the hood).*
 - Carefully wipe off the liquid around the sample and dropwise add 50 μ L antifade mounting solution to the samples. Add a coverslip to the samples using the blunt ends of tweezers and gently tap the coverslip to remove air bubbles.*
- 7. Observe and analyze with a fluorescence microscope.** CL488 is a green fluorescent dye with an excitation wavelength/emission wavelengths of 490/515 nm, respectively. CL594 is a red fluorescent dye with an excitation/emission wavelengths of 590/617 nm. Apoptotic cells are marked with bright green fluorescence. The negative control sample (those without TdT) should not have detectable signal.