

# Viability/Cytotoxicity cell staining kit (Calcein AM, EthD-1)

## Flow cytometry

1. Thaw reagents and prepare staining working solutions by adding 5  $\mu\text{L}$  of component A (Calcein AM) and 20  $\mu\text{L}$  of component B (EthD-1) to 10 mL of PBS, and vortex to mix.
2. Wash sample 2-3 times with 1X PBS through centrifugation.
3. Resuspend cells in 0.5 mL of staining working solution. The cell density should be  $1-5 \times 10^5$  cells/mL.  
*Note: We recommend preparing two tubes of additional samples, each tube with only one dye (Calcein AM and EthD-1) for compensation adjustment of each dye.*
4. Incubate for 15-20 min at room temperature, protected from light.
5. Assay the sample by flow cytometry. Calcein AM's excitation/emission maxima are 488/520nm and EthD-1's excitation/emission maxima are 528/617nm.  
*Note: Use single dye tubes to adjust and compensate, and the result should be two relatively independent cell populations: a live cell population with green fluorescence, and a dead cell population with red fluorescence.*

## Fluorescence microscopy

1. Thaw reagents and prepare staining working solutions by adding 5  $\mu\text{L}$  of component A (Calcein AM) and 20  $\mu\text{L}$  of component B (EthD-1) to 10 mL of PBS, and vortex to mix.  
*Note: For adherent cells, staining can be done directly. For suspension cells, centrifuge to collect cells for staining.*
2. Wash the cells 2-3 times with 1X PBS.
3. Aspirate and discard PBS. For adherent cells, add Calcein AM/ EthD-1 staining solution to cover the sample. For suspended cells, add enough working solution to bring cell density to  $1-5 \times 10^5$  cells/mL.
4. Incubate for 15-20 min at room temperature, protected from light.
5. Mount and observe the labeled cells under a fluorescence microscope.