



Ausgewählte Validierungsdaten



Immunofluorescent analysis of (-20°C Ethanol) fixed PC-3 cells using Coralite® Plus 488 S100A11 antibody (CL488-10237) at dilution of 1:200. Cell nuclei were stained with DAPI (blue).



1X10⁶ PC-3 cells were intracellularly stained with 0.4 ug Coralite® Plus 488 Anti-Human S100A11 (CL488-10237) (red), or 0.4 ug Control Antibody. Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).