

À des fins de recherche uniquement

# Anticorps Monoclonal anti-DEK

Numéro de catalogue: **CL647-66194** **Phare**



## Informations de base

Numéro de catalogue: CL647-66194	Numéro d'acquisition GenBank: BC035259	Méthode de purification: Purification par protéine A
Taille: 100ul , Concentration: 1000 µg/mL by Nanodrop;	Identification du gène (NCBI): 7913	CloneNo.: 5E3A6
Hôte: Mouse	Nom complet: DEK oncogene	Excitation/Emission maxima wavelengths: 654 nm / 674 nm
Isotype: IgG2a	MW calculé 43 kDa	
Immunogen Catalog Number: AG8225	MW observés: 40-50 kDa	

## Applications

Applications testées:  
FC (Intra)

Spécificité de l'espèce:  
Humain, rat, souris

## Informations générales

Dek oncogene(DEK) is first recognized in a fusion with CAN nucleoporin protein, and it's a also autoantigen in patients with pquciarticular onset juvenile rheumatoid arthritis. As a site-specific binding protein, DEK may involve in transcriptional regulation and signal transduction. Analysis by agarose gel electrophoresis, DEK showed to reduce the efficiency with which chromatin is replicated. It's one of the first components of the splicing complex, which remains associate with spliced exons dependent on prior splicing of pre-mRNA.

## Stockage

Stockage:  
Stocker à -20 °C. Éviter toute exposition à la lumière. Stable pendant un an après l'expédition.

Tampon de stockage:  
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.


L'aliquotage n'est pas nécessaire pour le stockage à -20C

**\*\*\* Les 20ul contiennent 0,1% de BSA.**

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com) W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



1X10<sup>6</sup> K-562 cells were intracellularly stained with 0.2 ug CoraLite® Plus 647 Anti-Human DEK (CL647-66194, Clone:5E3A6) (red), or 0.2 ug isotype control (blue). Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).