

À des fins de recherche uniquement

# Anticorps Recombinant de lapin anti- Phospho-RPS6KA1 (Ser380)

Numéro de catalogue:**CL647-80108**



## Informations de base

Numéro de catalogue:	BC014966	Méthode de purification:
CL647-80108		Purification par protéine A
Taille:	Identification du gène (NCBI):	CloneNo.:
100ul , Concentration: 1000 µg/ml by Nanodrop;	6195	7F23
Hôte:	Nom complet:	Excitation/Emission maxima
Lapin	ribosomal protein S6 kinase, 90kDa, polypeptide 1	wavelengths:
Isotype:	MW calculé	654 nm / 674 nm
IgG	735 aa, 83 kDa	
	MW observés:	
	83-90 kDa	

## Applications

Applications testées:

FC (Intra)

Spécificité de l'espèce:

Humain, souris

## Informations générales

The 90 kDa ribosomal S6 kinases (RSK1-4) are a family of widely expressed serine/threonine kinases and play an important role in the MAPK signaling cascade and act downstream of ERK1/2 signaling. It can mediate mitogenic and stress-induced activation of the transcription factors CREB1, ETV1/ER81 and NR4A1/NUR77. RSK1 (RPS6KA1) is also designated as p90RSK associated with cellular proliferation, survival and differentiation. It has been reported that the phosphorylation of RPS6KA1 at Ser380 is involved in autophosphorylation by activation of the C-terminal kinase domain. (PMID: 16458888, 30813401, 17259979, 10679322)

## Stockage

Stockage:

Stocker à -20 °C. Éviter toute exposition à la lumière.

Tampon de stockage:

PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.


L'aliquotage n'est pas nécessaire pour le stockage à -20°C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:  
T: 1(888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



1X10<sup>6</sup> Jurkat cells untreated (dashed lines) or treated with TPA (red) were intracellularly stained with 0.13 ug Coralite® Plus 647 Anti-Human Phospho-RPS6KA1 (Ser380) (CL647-80108, Clone:7F23) (red), or 0.13 ug Control Antibody (blue). Cells were fixed with 4% PFA and permeabilized with 90% MeOH.