# Recombinant GST Tag protein



# **Basic Information**

Catalog Number:

Ag25094

Size:

50 µg Form:

Available lyophilized

Species:

**Expression Source:** 

e coli.-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

MSPILGYWKIKGLVQPTRLLLEYLEEKYEEHLYERDEGD KWRNKKFELGLEFPNLPYYIDGDVKLTQSMAIIRYIADK HNMLGGCPKERAEISMLEGAVLDIRYGVSRIAYSKDFE TLKVDFLSKLPEMLKMFEDRLCHKTYLNGDHVTHPDF MLYDALDVVLYMDPMCLDAFPKLVCFKKRIEAIPQIDK YLKSSKYIAWPLQGWQATFGGGDHPPKSD

(1-220 aa encoded by FN315687)

# Reconstitution and Storage

#### Reconstitution:

Reconstitute at 0.25  $\mu g/\mu l$  in 200  $\mu l$  sterile water for short-term

After reconstitution with sterile water, if glycerol has no effect on subsequent experiments, it is recommended to add an equal volume of glycerol for long-term storage (see Stability and Storage for more details).

If a different concentration is needed for your purposes please adjust the reconstitution volume as required (please note: the ion concentration of the final solution will vary according to the volume used).

Note: Centrifuge vial before opening. When reconstituting, gently pipet and wash down the sides of the vial to ensure full recovery of the protein into solution.

### Shipping:

The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature (see below).

# **Purity**

85%, by SDS-PAGE with Coomassie Brilliant Blue staining.

# **Formulation**

The purified protein was Lyophilized from sterile PBS (58mM Na2HPO4,17mM NaH2PO4, 68mM NaCl, pH7.4). 5 %  $trehalose \ and \ 5\ \%\ mannitol\ are\ added\ as\ protectant\ before\ lyophilization.\ The\ elution\ buffer\ contain\ 300mM$ imidazole.

# Stability and Storage

Store for up to 12 months at -20°C to -80°C as lyophilized powder.

# Storage of **Reconstituted Protein**

**Short Term Storage:** Store at 2-8°C for (1-2 weeks).

Long Term Storage:

Aliquot and store at -20°C to -80°C for up to 3 months, reconstitution with sterile water and addition of an equal volume of glycerol. Avoid repeat freeze-thaw cycles.

# Selected Validation Data

