

#### For Research Use Only

# IHCeasy DDB1 Ready-To-Use IHC Kit

### Catalog Number: KHC1760

#### **General Information**

Sample type: FFPE tissue Cited sample type: Reactivity: Human, Mouse, Rat **Cited Reactivity:** 

Assay typ Immunohistochemistry Primary antibody type: Rabbit Polyclonal Secondary antibody type: Polymer-HRP-Goat anti-Rabbit

# **Kit Component**

Size	Concentration
100 mL	50×
100 mL ×2	20×
5 mL	RTU
5 mL	RTU
5 mL	RTU
0.2 mL	RTU
4 mL	RTU
5 mL	RTU
5 mL	RTU
5 mL	RTU
1 slide (Optional)	FFPE
1 Сору	
1 Сору	
	100 mL 100 mL ×2 5 mL 5 mL 5 mL 0.2 mL 4 mL 5 mL 5 mL 5 mL 5 mL 1 slide (Optional) 1 Copy

# **Storage Instructions**

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.

#### Background

DDB1, also named as XAP1, XPCe, DDBa and XPE-BF, belongs to the DDB1 family. It is required for DNA repair. DDB1 binds to DDB2 to form the UV-damaged DNA-binding protein complex (the UV-DDB complex). The UV-DDB complex may recognize UV-induced DNA damage and recruit proteins of the nucleotide excision repair pathway (the NER pathway) to initiate DNA repair. The functional specificity of the DCX E3 ubiquitin-protein ligase complex is determined by the variable substrate recognition component recruited by DDB1. The primary antibody in this kit is specific to DDB1.


#### **Synonyms**

DDB p127 subunit, DDB1, DDBA, DNA damage binding protein 1, DNA damage binding protein a, HBV X associated protein 1, UV damaged DNA binding factor, UV DDB 1, UV DDB1, XAP 1, XAP1, XPCE, XPE, XPE BF, XPE binding factor


For technical support and original validation data for this product please contact: T: 1 (888) 4PTGLAB (1-888-478-4522) (toll E: proteintech@ptglab.com free in USA), or 1(312) 455-8498 (outside W: ptglab.com USA)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffinembedded human stomach cancer tissue slide using KHC1760 (DDB1 IHC Kit).



Immunohistochemical analysis of paraffinembedded human appendicitis tissue slide using KHC1760 (DDB1 IHC Kit).



Immunohistochemical analysis of paraffinembedded mouse small intestine tissue slide using KHC1760 (DDB1 IHC Kit).




Immunohistochemical analysis of paraffinembedded mouse brain tissue slide using KHC1760 (DDB1 IHC Kit).



Immunohistochemical analysis of paraffinembedded mouse kidney tissue slide using KHC1760 (DDB1 IHC Kit).



Immunohistochemical analysis of paraffinembedded rat small intestine tissue slide using KHC1760 (DDB1 IHC Kit).



Immunohistochemical analysis of paraffinembedded rat kidney tissue slide using KHC1760 (DDB1 IHC Kit).