

For Research Use Only

IHCeasy TIPIN Ready-To-Use IHC Kit

Catalog Number: KHC2875

General Information

Sample type: FFPE tissue Cited sample type: Reactivity: Human, Mouse, Rat **Cited Reactivity:**

Assay typ Immunohistochemistry Primary antibody type: Rabbit Polyclonal Secondary antibody type: Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Сору	
Manual	1 Сору	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.

Background

TIPIN, also named as TIMELESS-interacting protein, is a 301 amino acid protein, which belongs to the CSM3 family. TIPIN Plays an important role in the control of DNA replication and the maintenance of replication fork stability. TIPIN forms a complex with TIMELESS and this complex regulates DNA replication processes under both normal and stress conditions, stabilizes replication forks and influences both CHEK1 phosphorylation and the intra-S phase checkpoint in response to genotoxic stress. stress.


Synonyms


TIMELESS interacting protein, TIMELESS-interacting protein


For technical support and original validation data for this product please contact: T: 1 (888) 4PTGLAB (1-888-478-4522) (toll E: proteintech@ptglab.com free in USA), or 1(312) 455-8498 (outside W: ptglab.com USA)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data








Immunohistochemical analysis of paraffinembedded human appendicitis tissue slide using KHC2875 (TIPIN IHC Kit).

Immunohistochemical analysis of paraffinembedded mouse small intestine tissue slide using KHC2875 (TIPIN IHC Kit).

Immunohistochemical analysis of paraffinembedded rat skin tissue slide using KHC2875 (TIPIN IHC Kit).



Immunohistochemical analysis of paraffinembedded rat small intestine tissue slide using KHC2875 (TIPIN IHC Kit).