



# IHCeasy® ZAK Ready-To-Use IHC Kit

Catalog Number: KHC2975

General Information

Sample type: FFPE tissue Cited sample type: Reactivity: Human Cited Reactivity: Assay type: Immunohistochemistry Primary antibody type: Rabbit Polyclonal

Secondary antibody type: Polymer-HRP-Goat anti-Rabbit

#### Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

### Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


### Background

ZAK(sterile-alpha motif and leucine zipper containing kinase AZK) is also named as MLTK, MAPKKK, mlklak, MLK7, AZK, MLT, MRK, HCCS-4, MAP3K20 and belongs to the MAPKKK family. It is a mitogenactivated protein kinase kinase kinase (MAP3K) that activates the stress-activated protein kinase/c-jun N-terminal kinase pathway and activates NF-kappaB. ZAK contributes to regulation of DNA damage checkpoints through a p38gamma-independent pathway.

## Synonyms

AZK, EC: 2.7.11.25, HCCS~4, HCCS-4, Leucine~zipper-~and~sterile~alpha~motif-containing~kinase

# Selected Validation Data



Immunohistochemical analysis of paraffinembedded human ovary cancer tissue slide using KHC2975 (ZAK IHC Kit).