

For Research Use Only

# STK16 Polyclonal antibody

Catalog Number: 10293-1-AP



## Basic Information

<b>Catalog Number:</b> 10293-1-AP	<b>GenBank Accession Number:</b> BC002618	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , Concentration: 350 µg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 8576	<b>Recommended Dilutions:</b> IP 0.5-4.0 ug for IP and 1:500-1:1000 for WB
<b>Source:</b> Rabbit	<b>Full Name:</b> serine/threonine kinase 16	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 39.4 kDa	
<b>Immunogen Catalog Number:</b> AG0339	<b>Observed MW:</b> 35 kDa	

## Applications

<b>Tested Applications:</b> IP, ELISA	<b>Positive Controls:</b> IP : mouse brain tissue,
<b>Species Specificity:</b> human, mouse	

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



IP result of anti-STK16(IP:10293-1-AP, 4ug;  
Detection:10293-1-AP 1:500) with mouse brain  
tissue lysate 1280 ug.