

For Research Use Only

TFDP1 Polyclonal antibody

Catalog Number: 11043-1-AP

1 Publications



Basic Information

Catalog Number:

11043-1-AP

Size:

150ul, Concentration: 400 ug/ml by Nanodrop and 253 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1511

GenBank Accession Number:

BC011685

GeneID (NCBI):

7027

UNIPROT ID:

Q14186

Full Name:

transcription factor Dp-1

Calculated MW:

46 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:20-1:200

Applications

Tested Applications:

IHC, ELISA

Cited Applications:

IF

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

IHC : human lymphoma tissue,

Note-IHC: suggested antigen retrieval with *TE buffer pH 9.0*; (*) Alternatively, antigen retrieval may be performed with *citrate buffer pH 6.0*

Notable Publications

Author	Pubmed ID	Journal	Application
Z L Liu	38519500	Nat Commun	IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human lymphoma using 11043-1-AP (TFDP1 antibody) at dilution of 1:50 (under 10x lens).