

For Research Use Only

UBE2K Polyclonal antibody

Catalog Number: 11834-3-AP

4 Publications



Basic Information

Catalog Number:

11834-3-AP

Size:

150ul, Concentration: 333 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2415

GenBank Accession Number:

BC022804

GeneID (NCBI):

3093

UNIPROT ID:

P61086

Full Name:

ubiquitin-conjugating enzyme E2K (UBC1 homolog, yeast)

Calculated MW:

200 aa, 22 kDa

Observed MW:

25 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

IHC 1:20-1:200

Applications

Tested Applications:

WB, IHC, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human, rat

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB: mouse brain tissue, mouse testis tissue, human brain tissue, rat brain tissue, rat testis tissue

IHC: human gliomas tissue, human thyroid cancer tissue

Notable Publications

Author	Pubmed ID	Journal	Application
Kazuhiko Watabe	32686212	Neuropathology	WB
Qingqing Wang	33647537	Biochem Biophys Res Commun	WB
Mo Juan Li	33685573	Biomed Environ Sci	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 11834-3-AP (UBE2K antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human gliomas using 11834-3-AP (UBE2K antibody) at dilution of 1:50 (under 10x lens).