

For Research Use Only

DDT Polyclonal antibody

Catalog Number: 12389-1-AP

1 Publications



Basic Information

Catalog Number:

12389-1-AP

Size:

150ul, Concentration: 800 ug/ml by Nanodrop and 400 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3065

GenBank Accession Number:

BC015508

GeneID (NCBI):

1652

UNIPROT ID:

P30046

Full Name:

D-dopachrome tautomerase

Calculated MW:

118 aa, 13 kDa

Observed MW:

13 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:300-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : rat liver tissue, HeLa cells

Notable Publications

Author	Pubmed ID	Journal	Application
Benjamin L Woolbright	36214539	J Pathol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat liver tissue were subjected to SDS PAGE followed by western blot with 12389-1-AP (DDT Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.