

For Research Use Only

SPOCK1 Polyclonal antibody

Catalog Number: 12512-1-AP

1 Publications



Basic Information

Catalog Number:

12512-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3129

GenBank Accession Number:

BC030691

GeneID (NCBI):

6695

UNIPROT ID:

Q08629

Full Name:

sparc/osteonectin, cwcw and kazal-like domains proteoglycan (testican)
1

Calculated MW:

439 aa, 49 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

IHC

Species Specificity:

human

Cited Species:

human

Notable Publications

Author	Pubmed ID	Journal	Application
Dehu Chen	28940639	J Cell Mol Med	IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

