

For Research Use Only

# SPOCK1 Polyclonal antibody

Catalog Number: 12512-1-AP **1 Publications**



## Basic Information

Catalog Number: 12512-1-AP	GenBank Accession Number: BC030691	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 6695	
Source: Rabbit	Full Name: sparc/osteonectin, cwcv and kazal- like domains proteoglycan (testican) 1	
Isotype: IgG	Calculated MW: 439 aa, 49 kDa	
Immunogen Catalog Number: AG3129		

## Applications

Tested Applications:  
ELISA

Cited Applications:  
IHC

Species Specificity:  
human

Cited Species:  
human

## Background Information

## Notable Publications

Author	Pubmed ID	Journal	Application
Dehu Chen	28940639	J Cell Mol Med	IHC

## Storage

Storage:  
Store at -20°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data