

For Research Use Only

MAPK6 Polyclonal antibody

Catalog Number: 12839-1-AP

Featured Product

3 Publications



Basic Information

Catalog Number:

12839-1-AP

Size:

150ul, Concentration: 400 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3577

GenBank Accession Number:

BC035492

GeneID (NCBI):

5597

UNIPROT ID:

Q16659

Full Name:

mitogen-activated protein kinase 6

Calculated MW:

83 kDa

Observed MW:

100-105 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human, rat

Positive Controls:

WB : A431 cells,

Notable Publications

Author	Pubmed ID	Journal	Application
Guanqun Huang	29764463	J Exp Clin Cancer Res	WB,IHC,IF
Chuanxian Hu	29456679	Exp Ther Med	WB
Juanjuan Tan	28429763	Sci Rep	WB,IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 12839-1-AP (MAPK6 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.