

For Research Use Only

# LILRB5 Polyclonal antibody

Catalog Number: 12881-1-AP



## Basic Information

<b>Catalog Number:</b> 12881-1-AP	<b>GenBank Accession Number:</b> BC025704	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , Concentration: 1000 µg/ml by Nanodrop and 433 µg/ml by Bradford method using BSA as the standard;	<b>GeneID (NCBI):</b> 10990	<b>Recommended Dilutions:</b> WB 1:500-1:1000
<b>Source:</b> Rabbit	<b>Full Name:</b> Leukocyte immunoglobulin-like receptor, subfamily B (with TM and ITIM domains), member 5	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 590 aa, 64 kDa	
<b>Immunogen Catalog Number:</b> AG3886	<b>Observed MW:</b> 64 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : COLO 320 cells, Jurkat cells
<b>Species Specificity:</b> human	

## Background Information

### Storage


**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12881-1-AP (LILRB5 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.