

For Research Use Only

# MSL3 Polyclonal antibody

Catalog Number: 12977-1-AP



## Basic Information

<b>Catalog Number:</b> 12977-1-AP	<b>GenBank Accession Number:</b> BC031210	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , Concentration: 1100 µg/ml by Nanodrop and 600 µg/ml by Bradford method using BSA as the standard;	<b>GeneID (NCBI):</b> 10943	<b>Recommended Dilutions:</b> WB 1:500-1:2000
<b>Source:</b> Rabbit	<b>Full Name:</b> male-specific lethal 3 homolog (Drosophila)	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 521 aa, 60 kDa	
<b>Immunogen Catalog Number:</b> AG3592	<b>Observed MW:</b> 60 kDa	

## Applications

**Tested Applications:**  
WB, ELISA

**Species Specificity:**  
human, mouse, rat

**Positive Controls:**

WB : mouse skeletal muscle tissue, human stomach tissue

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.


Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with 12977-1-AP (MSL3 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.