

For Research Use Only

FGF13 Polyclonal antibody

Catalog Number: 13201-1-AP

3 Publications



Basic Information

Catalog Number:

13201-1-AP

Size:

150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3942

GenBank Accession Number:

BC034340

GeneID (NCBI):

2258

UNIPROT ID:

Q92913

Full Name:

fibroblast growth factor 13

Calculated MW:

28 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Zhang Xiujun X	23554914	PLoS One	WB
Xiaojing Yue	29405966	Placenta	WB,IF
J-J Song	30720165	Eur Rev Med Pharmacol Sci	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

