

For Research Use Only

# CAMTA2 Polyclonal antibody

Catalog Number: 13470-1-AP



## Basic Information

|  |   |  |
|--|---|--|
| <b>Catalog Number:</b><br>13470-1-AP   | <b>GenBank Accession Number:</b><br>BC032381                      | <b>Purification Method:</b><br>Antigen affinity purification |
| <b>Size:</b><br>150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard; | <b>GeneID (NCBI):</b><br>23125                                    |  |
| <b>Source:</b><br>Rabbit   | <b>Full Name:</b><br>calmodulin binding transcription activator 2 |  |
| <b>Isotype:</b><br>IgG   | <b>Calculated MW:</b><br>1234 aa, 132 kDa                         |  |
| <b>Immunogen Catalog Number:</b><br>AG4263   |   |  |

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human, mouse, rat

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data