

For Research Use Only

ARHGAP20 Polyclonal antibody

Catalog Number: 13775-1-AP



Basic Information

Catalog Number:

13775-1-AP

Size:

150ul , Concentration: 187 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4764

GenBank Accession Number:

BC039340

GeneID (NCBI):

57569

UNIPROT ID:

Q9P2F6

Full Name:

Rho GTPase activating protein 20

Calculated MW:

133 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human brain tissue,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human brain tissue were subjected to SDS PAGE followed by western blot with 13775-1-AP (ARHGAP20 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.