

For Research Use Only

PKP4 Polyclonal antibody

Catalog Number: 14062-1-AP



Basic Information

Catalog Number:

14062-1-AP

Size:

150ul , Concentration: 300 ug/ml by Nanodrop and 227 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5113

GenBank Accession Number:

BC050308

GeneID (NCBI):

8502

UNIPROT ID:

Q99569

Full Name:

plakophilin 4

Calculated MW:

132 kDa

Observed MW:

125-140 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : A431 cells, HepG2 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 14062-1-AP (PKP4 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.