

For Research Use Only

IRGC Polyclonal antibody

Catalog Number: 14090-1-AP **1 Publications**



Basic Information

Catalog Number: 14090-1-AP	GenBank Accession Number: BC066939	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 240 µg/ml by Nanodrop and 160 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 56269	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: immunity-related GTPase family, cinema	
Isotype: IgG	Calculated MW: 50 kDa	
Immunogen Catalog Number: AG5279	Observed MW: 50 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse testis tissue, rat testis tissue
Cited Applications: WB	
Species Specificity: human, mouse, rat	
Cited Species: mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yuki Kaneda	35618043	Dev Biol	WB

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 14090-1-AP (IRGC antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.