

For Research Use Only

# VGLL4 Polyclonal antibody

Catalog Number: 15135-1-AP

1 Publications



## Basic Information

Catalog Number:

15135-1-AP

Size:

150ul, Concentration: 480 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7283

GenBank Accession Number:

BC000763

GeneID (NCBI):

9686

UNIPROT ID:

Q14135

Full Name:

vestigial like 4 (Drosophila)

Calculated MW:

31 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

## Background Information

VGLL4, also named as KIAA0121 and Vgl-4, is a co-transcriptional regulator with no previously known function in ES cells. It is a positive regulator for survival of human ES cells.

## Notable Publications

Author	Pubmed ID	Journal	Application
Qingzhe Wu	34597132	Sci Adv	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

