

For Research Use Only

VGLL4 Polyclonal antibody

Catalog Number: 15135-1-AP **1 Publications**



Basic Information

Catalog Number: 15135-1-AP	GenBank Accession Number: BC000763	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 480 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 9686	
Source: Rabbit	Full Name: vestigial like 4 (Drosophila)	
Isotype: IgG	Calculated MW: 31 kDa	
Immunogen Catalog Number: AG7283		

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human, mouse, rat

Cited Species:
human, mouse

Background Information

VGLL4, also named as KIAA0121 and Vgl-4, is a co-transcriptional regulator with no previously known function in ES cells. It is a positive regulator for survival of human ES cells.

Notable Publications

Author	Pubmed ID	Journal	Application
Qingzhe Wu	34597132	Sci Adv	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data