

For Research Use Only

FOSL2 Polyclonal antibody

Catalog Number: 15832-1-AP

Featured Product

8 Publications



Basic Information

Catalog Number:

15832-1-AP

Size:

150ul, Concentration: 1000 ug/ml by Nanodrop and 300 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8397

GenBank Accession Number:

BC008899

GeneID (NCBI):

2355

UNIPROT ID:

P15408

Full Name:

FOS-like antigen 2

Calculated MW:

326 aa, 35 kDa

Observed MW:

39 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

human, mouse, rat, pig

Positive Controls:

WB : MCF-7 cells,

Notable Publications

Author	Pubmed ID	Journal	Application
Yu Lian	30326930	J Exp Clin Cancer Res	WB
Xu Ye	29887867	Front Immunol	WB
Yuan Yuan Zhang	32377742	Mol Med Rep	WB, IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 15832-1-AP (FOSL2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.