

For Research Use Only

FOSL2 Polyclonal antibody

Catalog Number: 15832-1-AP

Featured Product

6 Publications



Basic Information

Catalog Number: 15832-1-AP	GenBank Accession Number: BC008899	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 300 µg/ml by Nanodrop and 300 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 2355	Recommended Dilutions: WB 1:200-1:1000
Source: Rabbit	Full Name: FOS-like antigen 2	
Isotype: IgG	Calculated MW: 326 aa, 35 kDa	
Immunogen Catalog Number: AG8397	Observed MW: 31 kDa	

Applications

Tested Applications: WB,ELISA	Positive Controls: WB : MCF-7 cells, Jurkat cells, PC-3 cells
Cited Applications: IHC, WB	
Species Specificity: human	
Cited Species: human, pig, rat	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yu Lian	30326930	J Exp Clin Cancer Res	WB
Xu Ye	29887867	Front Immunol	WB
Yuanyuan Zhang	32377742	Mol Med Rep	WB,IHC

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MCF-7 cells were subjected to SDS PAGE followed by western blot with 15832-1-AP (FOSL2 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.