

For Research Use Only

# FOXS1 Polyclonal antibody

Catalog Number: 16234-1-AP

1 Publications



## Basic Information

Catalog Number:

16234-1-AP

Size:

150ul, Concentration: 900 ug/ml by Nanodrop and 453 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8957

GenBank Accession Number:

BC013408

GeneID (NCBI):

2307

UNIPROT ID:

O43638

Full Name:

forkhead box S1

Calculated MW:

330 aa, 35 kDa

Observed MW:

35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : mouse liver tissue, rat liver tissue

## Notable Publications

Author	Pubmed ID	Journal	Application
Ruimin Ren	38780613	FASEB J	WB, IHC

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 16234-1-AP (FOXs1 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.