For Research Use Only

CENPL Polyclonal antibody

Catalog Number: 17007-1-AP

Featured Product

2 Publications



Purification Method:

WB 1:500-1:2400

Antigen affinity purification

Recommended Dilutions:

Basic Information

Catalog Number:

17007-1-AP

GenBank Accession Number:

BC033154

GeneID (NCBI):

Q8NoS6

Size: 150ul , Concentration: 133 ug/ml by 91687

Nanodrop and 133 ug/ml by Bradford UNIPROT ID:

method using BSA as the standard;

Source:

Full Name:

Rabbit centromere protein L Isotype:

Calculated MW: IgG 390 aa, 44 kDa

Immunogen Catalog Number: Observed MW:

AG10703 48 kDa

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls: WB: PC-3 cells, Raji cells

Notable Publications

Author	Pubmed ID	Journal	Application
Shan Wang	38522807	Cell Signal	WB
Kun He	37914022	Int J Biochem Cell Biol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

Selected Validation Data



PC-3 cells were subjected to SDS PAGE followed by western blot with 17007-1-AP (CENPL antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.