

For Research Use Only

SDR42E1 Polyclonal antibody

Catalog Number: 17347-1-AP



Basic Information

Catalog Number: 17347-1-AP	GenBank Accession Number: BC019670	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 300 µg/ml by Nanodrop and 207 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 93517	
Source: Rabbit	Full Name: short chain dehydrogenase/reductase family 42E, member 1	
Isotype: IgG	Calculated MW: 383 aa, 43 kDa	
Immunogen Catalog Number: AG11391		

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human bladder tissue slide using 17347-1-AP (SDR42E1 antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human bladder tissue slide using 17347-1-AP (SDR42E1 antibody) at dilution of 1:200 (under 40x lens).