

For Research Use Only

GPBP1L1 Polyclonal antibody

Catalog Number: 17375-1-AP



Basic Information

Catalog Number:

17375-1-AP

Size:

150ul , Concentration: 900 ug/ml by Nanodrop and 400 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11180

GenBank Accession Number:

BC022218

GeneID (NCBI):

60313

UNIPROT ID:

Q9HC44

Full Name:

GC-rich promoter binding protein 1-like 1

Calculated MW:

474 aa, 52 kDa

Observed MW:

58 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:5000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HepG2 cells, Jurkat cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17375-1-AP (GPBP1L1 antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours.