

For Research Use Only

WNT3 Polyclonal antibody

Catalog Number: 17983-1-AP

4 Publications



Basic Information

Catalog Number:

17983-1-AP

Size:

150ul , Concentration: 160 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12399

GenBank Accession Number:

BC114219

GeneID (NCBI):

7473

UNIPROT ID:

P56703

Full Name:

wingless-type MMTV integration site family, member 3

Calculated MW:

355 aa, 40 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human, mouse, pig

Notable Publications

| Author | Pubmed ID | Journal | Application |
|-------------|-----------|--------------|-------------|
| Qing Zuo | 27853262 | Sci Rep | WB |
| Jianqing Ye | 30431120 | Oncol Rep | WB |
| Min Li | 27247798 | J Ophthalmol | WB |

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

