

For Research Use Only

RHOBTB3 Polyclonal antibody

Catalog Number: 19014-1-AP

Featured Product

2 Publications



Basic Information

Catalog Number:

19014-1-AP

Size:

150ul, Concentration: 293 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5782

GenBank Accession Number:

BC041337

GeneID (NCBI):

22836

UNIPROT ID:

O94955

Full Name:

Rho-related BTB domain containing 3

Calculated MW:

611 aa, 69 kDa

Observed MW:

69 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Notable Publications

Author	Pubmed ID	Journal	Application
Maeve Long	32354068	Cells	WB
Lu Albert A	24145166	J Cell Biol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

