

For Research Use Only

TNFR2 / TNFRSF1B Polyclonal antibody

Catalog Number: 19272-1-AP

Featured Product

38 Publications



Basic Information

Catalog Number:

19272-1-AP

Size:

150ul, Concentration: 700 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5866

GenBank Accession Number:

BC052977

GeneID (NCBI):

7133

UNIPROT ID:

P20333

Full Name:

tumor necrosis factor receptor superfamily, member 1B

Calculated MW:

48 kDa

Observed MW:

70-75 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

IP 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

IHC 1:50-1:500

IF/ICC 1:50-1:500

Applications

Tested Applications:

WB, IHC, IF/ICC, IP, ELISA

Cited Applications:

WB, IHC, IF, IP, CoIP

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse, rat

Positive Controls:

WB : NK-92 cells, Jurkat cells, HEK-293 cells, THP-1 cells, mouse thymus tissue, rat thymus tissue

IP : HEK-293 cells,

IHC : human spleen tissue,

IF/ICC : THP-1 cells,

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Background Information

Tumor necrosis factor-alpha (TNFA/TNFSF2) is a multifunctional cytokine that plays a key role in regulating inflammation, immune functions, host defense, and apoptosis (PMID: 16407280). TNFA signals through two distinct cell surface receptors, TNFR1 (TNFRSF1A, CD120a, p55) and TNFR2 (TNFRSF1B, CD120b, p75). TNFR1 is widely expressed, whereas TNFR2 exhibits more restricted expression, being found on CD4 and CD8 T lymphocytes, endothelial cells, microglia, oligodendrocytes, neuron subtypes, cardiac myocytes, thymocytes and human mesenchymal stem cells (PMID: 20489699; 22374304). In contrast to TNFR1, TNFR2 does not have a death domain. TNFR2 only signals for antiapoptotic reactions. However, recent evidence indicates that TNFR2 also signals to induce TRAF2 degradation (PMID: 22374304). Various defects in the TNFR2 pathway, due to polymorphisms in the TNFR2 gene, upregulated expression of TNFR2 and TNFR2 shedding, have been implicated in the pathology of several autoimmune disorders (PMID: 20489699).

Notable Publications

Author	Pubmed ID	Journal	Application
Minami Uchida	31620105	Front Microbiol	IF
Di Huang	30224822	Nat Immunol	FC,IF
Qian Chen	30187338	Inflammation	

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 19272-1-AP (TNFR2 / TNFRSF1B antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.




Immunofluorescent analysis of (4% PFA) fixed THP-1 cells using TNFR2 / TNFRSF1B antibody (19272-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2).



IP result of anti-TNFR2 / TNFRSF1B (IP:19272-1-AP, 4ug; Detection:19272-1-AP 1:500) with HEK-293 cells lysate 1200ug.



Immunohistochemical analysis of paraffin-embedded human spleen tissue slide using 19272-1-AP (TNFR2 / TNFRSF1B Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human spleen tissue slide using 19272-1-AP (TNFR2 / TNFRSF1B Antibody) at dilution of 1:200 (under 40x lens).