

For Research Use Only

ENTPD4 Polyclonal antibody

Catalog Number: 19597-1-AP



Basic Information

Catalog Number:

19597-1-AP

Size:

150ul , Concentration: 300 ug/ml by Nanodrop and 287 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5175

GenBank Accession Number:

BC034477

GeneID (NCBI):

9583

UNIPROT ID:

Q9Y227

Full Name:

ectonucleoside triphosphate diphosphohydrolase 4

Calculated MW:

70 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : Jurkat cells, human placenta tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 19597-1-AP (ENTPD4 antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours.