

For Research Use Only

AFF1-Specific Polyclonal antibody

Catalog Number: 19975-1-AP



Basic Information

Catalog Number:

19975-1-AP

Size:

150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_005935

GeneID (NCBI):

4299

UNIPROT ID:

P51825

Full Name:

AF4/FMR2 family, member 1

Calculated MW:

131 kDa

Observed MW:

125-145 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

AFF1, also named as AF4, FEL, MLLT2 and PBM1, belongs to the AF4 family. A chromosomal aberration involving AFF1 is associated with acute leukemias. The antibody is specific to AFF1.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

