

For Research Use Only

ANKRD7 Polyclonal antibody

Catalog Number: 21023-1-AP



Basic Information

Catalog Number:

21023-1-AP

Size:

150ul, Concentration: 600 ug/ml by Nanodrop and 593 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15285

GenBank Accession Number:

BC032799

GeneID (NCBI):

56311

UNIPROT ID:

Q92527

Full Name:

ankyrin repeat domain 7

Calculated MW:

254 aa, 29 kDa

Observed MW:

31 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human testis tissue,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 21023-1-AP (ANKRD7 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.