

For Research Use Only

ESYT1 Polyclonal antibody

Catalog Number: 21450-1-AP



Basic Information

Catalog Number: 21450-1-AP	GenBank Accession Number: BC004998	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 300 µg/ml by Nanodrop;	GeneID (NCBI): 23344	Recommended Dilutions: WB 1:2000-1:8000
Source: Rabbit	Full Name: family with sequence similarity 62 (C2 domain containing), member A	
Isotype: IgG	Calculated MW: 1104 aa, 123 kDa	
Immunogen Catalog Number: AG13895	Observed MW: 123 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HEK-293T cells, HeLa cells, Jurkat cells
Species Specificity: human, mouse, rat	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 21450-1-AP (ESYT1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.