

For Research Use Only

TMTC1 Polyclonal antibody

Catalog Number: 23349-1-AP



Basic Information

Catalog Number: 23349-1-AP	GenBank Accession Number: BC042083	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 1000 µg/ml by Nanodrop and 520 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 83857	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: transmembrane and tetra-tricopeptide repeat containing 1	
Isotype: IgG	Calculated MW: 882 aa, 99 kDa	
Immunogen Catalog Number: AG19944	Observed MW: 99 kDa, 65 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse bladder tissue, mouse colon tissue, mouse spleen tissue
Species Specificity: human, mouse	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse bladder tissue were subjected to SDS PAGE followed by western blot with 23349-1-AP (TMTC1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.