

For Research Use Only

RBPJL Polyclonal antibody

Catalog Number: 24664-1-AP



Basic Information

Catalog Number:

24664-1-AP

Size:

150ul , Concentration: 650 ug/ml by Nanodrop and 387 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20284

GenBank Accession Number:

BC156875

GeneID (NCBI):

11317

UNIPROT ID:

Q9UBG7

Full Name:

recombination signal binding protein for immunoglobulin kappa J region-like

Calculated MW:

517 aa, 57 kDa

Observed MW:

57-69 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Jurkat cells, Raji cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 24664-1-AP (RBPJL Antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.