

For Research Use Only

RBPJL Polyclonal antibody

Catalog Number: 24664-1-AP



Basic Information

Catalog Number: 24664-1-AP	GenBank Accession Number: BC156875	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 650 µg/ml by Nanodrop and 387 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 11317	Recommended Dilutions: WB 1:500-1:3000
Source: Rabbit	Full Name: recombination signal binding protein for immunoglobulin kappa J region-like	
Isotype: IgG	Calculated MW: 517 aa, 57 kDa	
Immunogen Catalog Number: AG20284	Observed MW: 57-69 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Jurkat cells, Raji cells
Species Specificity: human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 24664-1-AP (RBPJL Antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.