

For Research Use Only

AWAT2 Polyclonal antibody

Catalog Number: 25119-1-AP



Basic Information

Catalog Number:

25119-1-AP

Size:

150ul , Concentration: 240 ug/ml by Nanodrop and 207 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18936

GenBank Accession Number:

BC153070

GeneID (NCBI):

158835

UNIPROT ID:

Q6E213

Full Name:

acyl-CoA wax alcohol acyltransferase

2

Calculated MW:

333 aa, 38 kDa

Observed MW:

38-42 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue, rat brain tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 25119-1-AP (AWAT2 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.