

For Research Use Only

SCML4 Polyclonal antibody

Catalog Number: 25439-1-AP



Basic Information

Catalog Number:

25439-1-AP

Size:

150ul , Concentration: 650 ug/ml by Nanodrop and 307 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22157

GenBank Accession Number:

BC033286

GeneID (NCBI):

256380

UNIPROT ID:

Q8N228

Full Name:

sex comb on midleg-like 4 (Drosophila)

Calculated MW:

414 aa, 45 kDa

Observed MW:

50-60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse liver tissue, A375 cells, HEK-293 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 25439-1-AP (SCML4 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.