

For Research Use Only

SCML4 Polyclonal antibody

Catalog Number: 25439-1-AP



Basic Information

| | | |
|--|---|--|
| Catalog Number: 25439-1-AP | GenBank Accession Number: BC033286 | Purification Method: Antigen affinity purification |
| Size: 150ul , Concentration: 650 µg/ml by Nanodrop and 307 µg/ml by Bradford method using BSA as the standard; | GeneID (NCBI): 256380 | Recommended Dilutions: WB 1:500-1:1000 |
| Source: Rabbit | Full Name: sex comb on midleg-like 4 (Drosophila) | |
| Isotype: IgG | Calculated MW: 414 aa, 45 kDa | |
| Immunogen Catalog Number: AG22157 | Observed MW: 50-60 kDa | |

Applications

| | |
|---|---|
| Tested Applications: WB, ELISA | Positive Controls: WB : mouse liver tissue, A375 cells, HEK-293 cells |
| Species Specificity: human, mouse | |

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 25439-1-AP (SCML4 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.